

▪ Raise the Standard ▪ Cystoscopy - Hysteroscopy



Performance

Our rigid scopes offer superior imaging. Higher light fiber density ensures more light at the operative site, while our sophisticated lens systems provide optimum performance, image quality and optical efficiency.

Bring out the Best of your HD Camera

Paired with your HD Camera systems, our scopes provide superior visualization in the OR.

Durability

Featuring soldered Sapphire scratch resistant windows and innovative rigid rod lens systems, ensuring the reliability you need.

Autoclavability

MediVision endoscopes are Guaranteed for autoclavability for 1 year.

Compatibility

Seamlessly compatible with Storz® and Linvatec® style cannula systems.

Reliability

All MediVision endoscopes are designed and manufactured in the USA..

Product Ordering Information

Cystoscopes / Hysteroscopes

MV400	4.0 mm, 0°, 302 mm working length
MV412	4.0 mm, 12°, 302 mm working length
MV430	4.0 mm, 30°, 302 mm working length
MV470	4.0 mm, 70°, 302 mm working length
MV2900	2.9 mm, 0°, 302 mm working length
MV2930	2.9 mm, 30°, 302 mm working length

Feature:	Benefit:
Soldered eyepiece window	Eliminates leaks, providing a longer working life for the scope
Thermal expansion compensator	Assures repeatable autoclavability
Sapphire objective window and eyepiece	Eliminates scratching
Maximum fiber packing capacity	Provides greater light at operative site - More fibers in scope
Color coded EZ-ID™ eyepieces	Provides easy identification - No need for ID tape
US manufacturer to FDA quality system	Guarantees quality

Service

Service as well as repair exchange programs available for all MediVision scopes.

MediVision also offers a high quality, economical and rapid solution to your endoscopic service needs. Our technicians service over 750 different models of rigid endoscopes from all major manufacturers.

**4883 E. La Palma Ave.,
Suite 503
Anaheim, CA 92807
Phone (714) 563-2772
Fax (714) 563 2711
www.medivisionusa.com**

The information presented in this brochure is intended to demonstrate a MediVision product. Always refer to the package insert, product label and/or user instructions before using any MediVision product. Products may not be available in all markets. Product availability is subject to the regulatory and medical practices that govern individual markets. Please contact your MediVision representative if you have questions about the availability of MediVision products in your area.

Products referenced with ™ designation are trademarks of MediVision.

Products referenced with ® designation are registered trademarks of MediVision. Storz® is a registered trademark of Karl Storz and Linvatec® is a registered trademark of Linvatec.